DoseEdge

October, 2012

What is it? Limitations?

- A system to improve patient safety in preparing IV solutions, reducing waste, and improved nurse tracking (minimizing telephone calls to get IV's)
- Right now, this system is unavailable to float nurses, or nurses not permanently assigned to a unit due to HIPAA concerns
- Will not be available on units that are not a location in Cerner, ex OR, PACU, Dialysis, Cath lab

Why?

- This is a means to provide communication and transparency to nursing as to the process of IV medications being prepared.
- This is also a means to reduce the number of telephone calls from nursing.
- The times on the screen are the administration times identified in Cerner, so if it requires changing, send a MedRequest.

Accessing DoseEdge

From your computer desktop, double click the icon displayed below




Enter e-mail User Name and Password




DoseEdge® is not intended to replace the knowledge, judgment or expertise of pharmacists and pharmacy technicians in the preparation of IV admixtures and oral doses. DoseEdge is a software system intended to improve and document workflow in your pharmacy, and help in reducing opportunities for errors. While the software can provide guidance that assists the pharmacist or technician in proper preparation of IV admixtures and oral doses, it cannot, and does not, make any decision regarding dose preparation or assure that the personnel use correct technique or prepare doses in the proper order. Pharmacists and technicians are expected to use DoseEdge in combination with proper training, supplies and techniques and to assume responsibility for the doses they prepare. Baxa Corporation and its suppliers disclaim all liability arising from the failure of the pharmacist or technician to use the proper standard of care in the preparation and handling of IV admixtures and oral doses. Use of the DoseEdge service and all related software is subject to all of the terms and conditions of your institution's DoseEdge Service Agreement. Patent: 7,096,212


Select "View Status Board"




The Status Board

Various stages of IV preparation








Box 3 DOSE DGE **Aultman Health Foundation** Freeze Updates MA5 Status Board MA5 Status Board Dose Legend:DC, Stat, First Dose, Dilution, Stock, Batch (Minutes indicate time the dose is in this MedRequest 1 hr before needed (Total: 1) Doses Waiting to be Prepared (Total: 3) Doses in Process (Total: 0) 640 Minutes 42 Minutes hydromorphone 10mg (MA5 - 5556-A), PatientID: acetaminophen 1,000mg (MA5 - 5524-A), PatientID Admin:00:00 42 Minutes acetaminophen 1,000mg (MA5 - 5504-A), PatientID: Admin:00:00 42 Minutes acetaminophen 1,000mg (MA5 - 5566-A), PatientID: Admin:00:00 Doses Waiting for Final Label (Total: 3) Doses Ready for Delivery (Total: 19) Doses Delivered to MA5 (Total: 10) ADK 4 Minutes BUR 15 Minutes 324 Minutes Admin:23:00 ampicillin-Admin:10:00 ceFAZolin 1g (MA5 - 5570-A), PatientID: The upper right box indicates items currently being A), Patien 18 Minutes 123 Minutes prepared. The time indicates how long those orders ceFAZolin 2g (MA5 - 5510-A), PatientID: ide 0.9% Soln 50mL (MA5 - 5512have been in that status as well. 17 Minutes 119 Minutes 261 Minutes Admin:00:00 ceFAZolin 1g (MA5 - 5552-A), PatientID: ampicillin-sulbactam 3g,0.9% sodium chloride 100ml in 100mL (MA5 - 5530dextrose 5% in 0.45% nacl 1,000mL, potassium chloride 20mEq (MA5 - 5512-A), PatientID: Admin:18:00 112 Minutes 240 Minutes clindamycin 600mg, Dextrose 5% in Water Addva 50mL (MA5 - 5502-A), vancomycin 1.25g, Sodium Chloride 0.9% 250mL (MA5 - 5518-A), PatientID: Admin:18:00 Admin:12:00 144 Minutes ceFAZolin 2g (MA5 - 5510-A), PatientID: Admin:16:30 Lactated Ringers 1,000mL, potassium chloride 20mEq (MA5 - 5528-A).

What's going on?

The technician preparing the IV is required to scan each component and take a picture of the intermediate steps.







Picture taking

Final picture – if you receive a product with NOT CHECKED, return to Pharmacy





2 barcodes on final label

Cerner scan




Pharmacy scan






ScanEvents

Tubed (with time tubed)



Or Delivered (with time delivered)



Questions

- Refer any questions to:
 - Jim Toohey, M.S., R.Ph.
 - jtoohey@aultman.com
 - 3-2877