


# ELECTROLYSIS OF WATER LAB SET UP

## EQUIPMENT



## J- HOOK PLACEMENT



J HOOKS ARE PLACED ON THE **OUTSIDE** OF THE TEST TUBE, AND IN THE NOTCH IN THE TEST TUBE HOLDERS, AS SHOWN ABOVE LEFT AND RIGHT.

BE SURE THAT J HOOK TERMINAL IS UNDERNEATH AND BELOW THE MOUTH OF THE TEST TUBE, AS SHOWN TO THE LEFT, **NOT INSIDE THE TEST TUBE!!!**