

Student: _____

Date: _____

Steps 35-1 Preparing a Syringe for Use

	S	U
1. Carries out Standard Steps A, B, C, D, and E as need indicates.	<input type="checkbox"/>	<input type="checkbox"/>
2. Removes outer wrapper or casing or assembles cartridge holder.	<input type="checkbox"/>	<input type="checkbox"/>
3. Removes excess air from cartridge if applicable.	<input type="checkbox"/>	<input type="checkbox"/>
4. Attaches a sterile needle to the syringe if needed without contaminating either item.	<input type="checkbox"/>	<input type="checkbox"/>
5. Maintains sterility of syringe.	<input type="checkbox"/>	<input type="checkbox"/>
6. Carries out Standard Steps X, Y, and Z.	<input type="checkbox"/>	<input type="checkbox"/>

Successfully completed: ☐

Needs practice and retesting: ☐

Comments:

Instructor: _____