

Steps 22-1 Weighing the Adult with a Standing Balance Scale

Student: _____

Date: _____

	S	U
1. Checks that scale calibration is at "0."	<input type="checkbox"/>	<input type="checkbox"/>
2. Places cover on footplate of scale.	<input type="checkbox"/>	<input type="checkbox"/>
3. Assists patient onto the scale.	<input type="checkbox"/>	<input type="checkbox"/>
4. Sets large weight indicator.	<input type="checkbox"/>	<input type="checkbox"/>
5. Adjusts weights until beam balances.	<input type="checkbox"/>	<input type="checkbox"/>
6. Records the weight.	<input type="checkbox"/>	<input type="checkbox"/>
7. Assists patient off of the scale.	<input type="checkbox"/>	<input type="checkbox"/>

Successfully completed ☐

Needs practice and retesting ☐

Comments:

Instructor: _____